

Instructions For Use

De Mayo D²® Knee Positioner®

Catalog Numbers	803-114, 803-114S		Status	APPROVED
Effective Date	May 26, 2020	Version	3	Page 1 of 2

1. Replace OR Table Pads with green *imp*® Pads.

Note: It is not necessary to remove X-Ray cassette plates.



2. Position the patient with the gluteal fold at the separation of the pads.



3. During Prep of the surgical leg, remove single OR Table Pad.



4. Create a 3"-5" fold in the final drape under the buttocks.



5. Place the De Mayo D²® Knee Positioner® in the well to achieve maximum flexion.



6. Check flexion for final approval of freedom of the drape to allow full range of motion, when the positioner is locked on the drapes.



7. Position clamp directly under the guide pins in accordance to the clamp IFU

See Individual Clamp IFUs

8. Insert *imp*[®] Patient Protective Pad[®] into the sterile boot

For patient's safety, always use *imp*[®] Patient Protective Pads[®].



9. Wrap cohesive bandage around the foot starting with a minimum of six (6) foot wraps, tear, and finish the wrap above the distractor block around the calf.

Note: Do not wrap cohesive bandage over the distractor block on the back of the boot.




Safety Test



Tighten clockwise to lock the knob and lever

1. Check all knobs and levers are locked
2. Patient is fully protected by *imp*[®] foam and cohesive wrap
3. Patient's leg can reach full flexion
4. Check that both Guide Pins are fully seated



Scan for additional documentation



De Mayo D2[®] Knee Positioner is a registered trademark of Innovative Medical Products, Inc.
imp[®] products are protected by patent & patent pending rights ~ go to www.imppatents.com
 © May 2020 Innovative Medical Products, Inc.
 ALL RIGHTS RESERVED



Innovative Medical Products
 87 Spring Lane
 Plainville, CT 06062
 P: 860-793-0391
 F: 866-459-1805