

Info Sheet

Knee Positioner Serialization Locations

Catalog Numbers	803, 803-S, 803-114, 803-114-S, 114, 114-S	Status	APPROVED
Effective Date	May 7, 2020	Version	2
		Page 1 of 1	

Keep all Critical Parts Together

Laser etched numbers assure that the correct components (*base, carriage, boot and handle*) stay together. Extends life, reduces wear and simplifies assembly.


All critical parts are Serialized with the same unique 3-digit code

1. Carriage
2. Boot
3. Base plate
4. Handle (D2 and V2 Only)

D2® (CAT# 803-114)

V2® (CAT# 114)

De Mayo (CAT# 803)



Scan for additional documentation

