

imp Sterilization Case Systems


Designed to accommodate the
De Mayo Knee Positioners® & Accessories

● 706-S will accept *any* Baseplate or the V² E™ Extension Arm



De Mayo Knee Positioner®


706-S
32" x 12" x 2.5"



D²




V²




V² E™ ~Extension

● Helps protect equipment




● Low profile cases comply with hospital weight guidelines

● Improves handling efficiency in the OR's sterile field



706-XS

14.5" x 12" x 6.5"



706-M

16.5" x 12" 6.5"

Various Kit Options Available!

imp® Sterilization Case Systems

Ordering Information

<i>Description</i>	<i>Catalog No.</i>
imp® X-Small Sterilization Case for Knee Positioner® Boot & Clamp Dimensions: 14.5" x 12" x 6.5" Weight: 10 lbs w/equipment	706-XS
imp® Small Sterilization Case for <i>any</i> De Mayo Knee Positioner® Baseplates Dimensions: 32" x 12" x 2.5" Weight: 13 lbs Base only ~ or Extension Arm	706-S
imp® Medium Sterilization Case for Knee Positioners® Boot, Clamp & Distractor Dimensions: 16.5" x 12" x 6.5" Weight: 15 lbs w/equipment	706-M
imp® Sterilization Case Kit for De Mayo Knee Positioners® Includes: 1 each ~ 706-XS and 706-S	706-XS KIT
imp® Sterilization Case Kit for De Mayo V ² E™ Knee Positioner® Includes: 2 ~ 706-S and 706-XS	706-XS-EXT KIT
imp® Sterilization Case Kit for De Mayo Knee Positioners® Includes: 1 each~ 706-M and 706-S	706-M KIT
imp® Sterilization Case Kit for De Mayo V ² E™ Positioners® Includes: 2~ 706-S and 1~ 706-M	706-M-EXT KIT
imp® Promo Sterilization Case for De Mayo Knee Positioners® Includes: 706-XS and 706-S, 1210 ASSEMBLY KIT and 1210-SP	706-XS-PROMO KIT
imp® Promo Sterilization Case for De Mayo Knee Positioners® Includes: 706-S, 706-M, 1210 ASSEMBLY KIT and 1210-SP	706-M-PROMO KIT



**Positioners & Additional Accessories
Sold Separately**



De Mayo Knee Positioner® is a registered trademark of Innovative Medical Products, Inc.
U.S. Patent No. 7,380,299
De Mayo V²E™ is a trademark of Innovative Medical Products, Inc.
U.S. Patent No. 8,132,278
De Mayo Universal Distractor® is a registered trademark of Innovative Medical Products, Inc.
U.S. Patent No. 8,048,082
© September 2017
ALL RIGHTS RESERVED

***imp*® Innovative**
Medical Products, Inc.
87 Spring Lane, PO Box 8028
Plainville, CT 06062
PH: (860) 793-0391 FAX: (866) 459-1805
Toll Free PH: 1-800-467-4944
www.impmedical.com