


Info Sheet

De Mayo Adapt2Fit™ Modular Knee Positioner Sterilization Tray Configuration

Catalog Number	1020	Effective Date	November 12, 2021	Status	APPROVED
Doc ID	INFO-IMP-0006	Version	2	Page 1 of 1	



1. Right baseplate
2. Left baseplate
3. Carriage
4. Baseplate connector
5. Boot
6. Quad[®] Clamp
7. Guide pins



Scan for additional documentation



imp[®] products are protected by patent & patent pending rights ~ go to www.imppatents.com
 © November 2021 Innovative Medical Products, Inc.
 ALL RIGHTS RESERVED



Innovative Medical Products
 87 Spring Lane
 Plainville, CT 06062
 P: 860-793-0391
 F: 866-459-1805